

The MOST ADVANCED Medical Alert System On The Market Today



What is a Personal Emergency Response System (PERS)?

MediPendant's™ personal emergency response system (PERS), is the medical alert device that provides professional assistance at the touch of a button. There are no telephone numbers to look up, and no need to dial the phone. With MediPendant's™ two-way voice technology, simply press the button, and communicate with an EMT-certified operator directly through the pendant.

How does it work?

When the large button on the pendant is pressed, a telephone call is made to an EMT-certified emergency operator.

Then what happens?

The operator will immediately dispatch medical personnel or a loved one to your location. They will also contact anyone that you have listed on your custom call list.



Be Assured


If you are concerned about a loved one, a friend, or anyone who is living alone, or may need assistance in an emergency, then it's time for MediPendant™. It will provide you with the most important thing... *peace of mind.*

COMPARISON CHART

Medical Alarm Features	MediPendant™	Competitor*
Speak and Listen Directly Through the Pendant	YES	No
Up to 600 + feet of Communication Range	YES	No
Out of Range Notification	YES	No
Call Center Staffed with Certified EMTs	YES	No
Call Center Guaranteed To Stay On Line Until Help Arrives	YES	No
Wallet Card ID Program	YES	No
Spanish Speaking Equipment	YES	No
Panic Button Range	Up to 600+ ft	300 ft
Backup Battery Life	Up to 30 hrs	22 hrs

* based on latest competitor equipment and pricing

In an emergency, time is of the utmost importance. If you ever require assistance, the MediPendant™ personal emergency response system will summon the help you need.



For more information, contact your local dealer.

Features

- Speak and listen to the operator directly through the pendant
- Easy installation/User friendly
- Operates anywhere from 300 feet to 600 feet* in or around the home †
- Water-resistant – may be worn in bath or shower**
- Can be worn around the neck, clipped to your belt, or around your wrist.
- Battery provides up to six hours of talk time and up to one year of stand-by time.
- 24/7 monitoring with EMT-certified professionals
- Bi-lingual equipment and operators available

*Line of sight.

† Building structure and conditions may affect range.

**Base station should never be exposed to water.

Your MediPendant™ System Includes:

- Two-way voice speakerphone pendant
- Easy to follow set-up guide
- Belt, neck and wrist attachments
- 1 battery backup
- 1 phone connection kit
- 1 power supply



Questions & Answers

Q: What is the advantage of MediPendant™?

A: MediPendant™ enables you to **speak and listen directly through the pendant.** There is no need to be near a base station.

Q: Is MediPendant™ difficult to install?

A: It is very simple to install. Just plug in the power cord and the phone jack. There are no tools required.

Q: Will my phone bill increase?

A: No, the system dials a pre-programmed toll-free number.

Q: What if there is a loss of power?

A: The base station has backup battery power which allows the system to be used if the power fails.

Q: Who answers the call when the button is pressed?

A: An emergency response call center operator will answer the call. All of our operators are EMT-certified and trained to provide on-call medical assistance.

Q: What happens if I am unable to speak?

A: When the button is pressed, you will hear a voice recording letting you know that your call for help has gone through. Within seconds, an operator will come on the line. If the operator cannot hear you, help will be dispatched immediately.


HOW DO I ORDER?

Contact your nationwide dealer.

DigiKnox
1831 Truxtun Ave. Suite 115
Bakersfield, CA 93301

877-322-4204
medipendantmedicalalarm.com

**The MEDICAL ALERT SYSTEM
that allows you to
SPEAK and LISTEN DIRECTLY
THROUGH THE PENDANT**



MediPendant™
Experience the Difference